

Brussel, 7 April 2025

Kwaliteit van Laboratoria

GLOBAL RAPPORT EKE HbA1c VERS BLOED - 2025.1C

Target HbA1c = 57.4 mmol/mol (Referentiemethode NGSP/DCCT)

INSTRUMENT	RESULT (Mean)	N	SD	CV%	Bias	Bias %
ADAMS HA-8180T IFCC mmol/mol	59.5	2	2.12	3.6	2.1	3.7
ADAMS HA-8180V IFCC mmol/mol	57.5	21	1.26	2.2	0.1	0.2
ADAMS HA-8190V IFCC mmol/mol	59.0	2	1.41	2.4	1.6	2.8
Afinion 2 IFCC mmol/mol	57.0	1	****	****	-0.4	-0.7
Capillarys 2 flex-piercing IFCC mmol/mol	57.0	1	****	****	-0.4	-0.7
Capillarys 3 Octa IFCC mmol/mol	56.3	11	1.27	2.3	-1.1	-1.9
Capillarys 3 Tera IFCC mmol/mol	56.4	5	1.14	2.0	-1.0	-1.7
cobas b 101 IFCC mmol/mol	55.0	1	****	****	-2.4	-4.1
cobas c 501/502 cobas 6000 8000 IFCC	55.0	2	4.17	7.6	-2.4	-4.2
D10 A1c IFCC mmol/mol	58.7	2	1.90	3.3	1.3	2.2
D100 A1c IFCC mmol/mol	57.9	8	1.38	2.4	0.5	0.9
Hb NEXT IFCC mmol/mol	61.4	4	1.70	2.8	4.0	7.0
Minicap flex-piercing IFCC mmol/mol	55.4	5	1.51	2.7	-2.0	-3.4
TOSOH G11 IFCC mmol/mol	58.8	28	0.70	1.2	1.4	2.5
TOSOH G8 IFCC mmol/mol	58.7	11	0.80	1.4	1.4	2.4

Uw individueel rapport kunt u opvragen via de site
https://qseenet.com/CueSee_frame.asp?SelectMenu1=2

De enquêtecoördinatoren
Y. Lenga

Dr. K. Vernelen/ A. Vantorre